

Device compatibility



Check if your smartphone is compatible with your Rexton devices

Direct streaming to hearing aids via Bluetooth

Apple devices:


Rexton Mfi (made for iPhone, iPad or iPod touch) hearing aids connect directly to your iPhone, iPad or iPod so you can stream your phone calls and music directly into your hearing aids.

Android devices:

With Rexton BiCore devices, you can now also stream directly to Android devices via the ASHA (Audio Streaming for Hearing Aids) standard.

ASHA-supported devices:

- Samsung Galaxy S21
- Samsung Galaxy S21 5G (SM-G991U)(US)
- Samsung Galaxy S21 (US)
- Samsung Galaxy S21+ 5G (SM-G996U)(US)
- Samsung Galaxy S21 Ultra 5G (SM-G998U)(US)
- Samsung Galaxy S21 5G (SM-G991B)
- Samsung Galaxy S21+ 5G (SM-G996B)
- Samsung Galaxy S21 Ultra 5G (SM-G998B)
- Samsung Galaxy Note 20 Ultra (SM-G)
- Samsung Galaxy Note 20 Ultra (SM-G)(US)
- Samsung Galaxy S20+ (SM-G)
- Samsung Galaxy S20+ (SM-G) (US)
- Samsung Galaxy S20 5G (SM-G981B)
- Samsung Galaxy S20 5G (SM-G981U1) (US)
- Samsung Galaxy S20 Ultra 5G (SM-G988B)
- Samsung Galaxy S20 Ultra 5G (SM-G988U)(US)
- Samsung Galaxy S20 (SM-G980F)
- Samsung Galaxy S20 (SM-G) (US)
- Samsung Galaxy Note20 5G (SM-N981U1) (US)
- Samsung Galaxy Note 10+ (SM-N975F)
- Samsung Galaxy Note 10+ (SM-N975U1)(US)



- Samsung Galaxy Note 10 (SM-N970F)
- Samsung Galaxy Note 10 (SM-N970U)(US)
- Samsung Galaxy Note 10 Lite (SM-N770F/DS)
- Samsung Galaxy S10 Lite (SM-G770F/DS)
- Samsung Galaxy S10 (SM-G973F)
- Samsung Galaxy S10 (SM-G973U1) (US)
- Samsung Galaxy S10+ (SM-G975F)
- Samsung Galaxy S10+ (SM-G975U) (US)
- Samsung Galaxy S10e (SM-G970F)
- Samsung Galaxy S10e (SM-G970U) (US/China)
- Samsung Galaxy S10 5G (SM-G977U) (US)
- Samsung Galaxy S10 5G (SM-G977B)
- Samsung Galaxy S9 (SM-G960F)
- Samsung Galaxy S9 (SM-G960U1) (US)
- Google Pixel 5
- Google Pixel 4a 5G
- Google Pixel 4
- Google Pixel 4a
- Google Pixel 4 XL
- Google Pixel 3XL
- Google Pixel 3
- Google Pixel 3a
- Google Pixel 3a XL
- OnePlus Nord

Mobile devices compatible with the Rexton app:



- **Apple devices:** The Rexton app is compatible with all Apple devices with the iOS version 11.0 and higher.
- **Android devices:** The Rexton app is compatible with most of the Android (Android 8.0 and higher) devices.

For the optimal stability of the app while using with our Bluetooth®, hearing aids we recommend using the following Android devices:

Samsung Galaxy S10+ (Android 9.0 or higher)

- SM-G975F

Samsung Galaxy S10 (Android 9.0 or higher)

- SM-G973F

Samsung Galaxy S9 (Android 8.0 or higher)

- SM-G960F
- SM-G960U1

Samsung Galaxy S8 (Android 7.0 or higher)

- SM-G950F
- SM-G950U1

Samsung Galaxy Note 8 (Android 7.1.1 or higher)

- SM-N950U1

Samsung Galaxy Note 4 (Android 6.0.1 or higher)

- SM-N910F

Samsung Galaxy S7 (Android 7.0 or higher)

- SM-G930FD
- SM-G930V

Samsung Galaxy S7 Edge (Android 7.0 or higher)

- SM-G935F

Samsung Galaxy J7 (Android 7.0 or higher)

- SM-J730F

Samsung Galaxy S6 (Android 6.0.1 or higher)

- SM-G920I

Samsung Galaxy S6 Edge (Android 7.0 or higher)

- SM-G925F

Samsung Galaxy S5 (Android 6.0.1 or higher)

- SM-G900F
- SM-G900V

Samsung Galaxy A5 (Android 8.0 or higher)

- SM-A520F
- SM-A510F

Samsung Galaxy A3 (Android 8.0 or higher)

- SM-A320FL

Samsung Galaxy A80 (Android 9.0 or higher)

- SM-A805F

Samsung Galaxy A50 (Android 9.0 or higher)

- SM-A505F

Samsung Galaxy A40 (Android 9.0 or higher)

- SM-A405FN

Google Pixel 3 (Android 9.0 or higher)

- All models

Google Pixel 2 (Android 9.0 or higher)

- All models

Google Pixel (Android 7.1.2 or higher)

- All models

Google Pixel XL (Android 9.0 or higher)

- All models

Huawei P20 (Android 8.1 or higher)

- EML-29

Huawei Nexus 6P (Android 8.0 or higher)

- Nexus 6P

OnePlus 7 (Android 9.0 or higher)

- GM1903

OnePlus 6 (Android 8.1.0 or higher)

- A6003

Nokia 8.1 (Android 9 or higher)

- TA-1121